AirSeal[®] i.F.S.

Intelligent Flow System


Instructions for use **EN**

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EN

This manual contains information that is subject to copyright. All rights reserved. This manual should not be photocopied, duplicated on microfilm or otherwise copied or distributed, completely or in part, without the approval of SURGIQUEST.

Some of the parts and equipment referred to in this manual bear registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark indicates that any given designation is not subject to trademark protection.

Users of SURGIQUEST products should not hesitate to point out to us any errors or unclari- ties in this manual.

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Table of contents

1	Important User Notes	3
2	Safety Information.....	4
2.1	Hazards.....	4
3	Device Purpose.....	7
3.1	Device-inherent Dangers	7
4	Initial Device Startup	11
4.1	Gas Connection	11
4.1.1	Connecting a Gas Cylinder.....	12
4.1.2	Connecting to Central Gas Supply.....	12
4.1.3	Gas Consumption Display.....	12
5	Operating the Device	14
5.1	Device Front	14
5.2	Device Rear	14
5.3	Bottom of the Device	14
5.4	Display	15
5.5	Switch Device On	15
5.5.1	Selecting Operating Mode.....	16
5.5.2	Insufflation Tube Sets	16
5.5.3	Starting/Stopping Insufflation.....	17
5.5.4	Power Device Off	18
6	Using and Controlling the AirSeal® i.F.S. in its Different Modes	19
6.1	Standard Insufflation Mode	19
6.2	Smoke Evacuation Mode.....	21
6.3	AirSeal Mode.....	22
6.3.1	AirSeal® Optical Obturator and Cannula System.....	23
6.3.2	Initial Insufflation.....	24
7	Safety Functions.....	27
7.1	General Safety Functions	27
7.2	Contamination Alarm	27
7.3	Safety Functions in AirSeal Mode	28
7.4	Fill Level Display	28
8	User Menu	29
8.1	Setting First Nominal Pressure.....	30
8.2	Gas Flow Rates	30
8.3	Setting Smoke Evacuation Level	30
8.4	Setting General Volume	30
8.5	Setting Brightness.....	30
8.6	Setting Audible Occlusion Signal.....	31
8.7	Setting Gas Supply	31
8.8	Setting the Language.....	31
8.9	Checking Software Version.....	31
8.10	Resetting or Restoring to Factory Setting	31
9	Care and Maintenance	32
9.1	Cleaning the Device.....	32
9.2	Annual Inspection	32
9.3	Maintenance by Authorized Service Technician.....	32
9.4	Replacing the Fuse	33
10	Annual Inspection	34
10.1	Electrical Safety Test	34
10.2	Basic Function Test.....	34
10.3	Pressure Sensor Test	35
10.4	Pressure Monitoring Test	36
10.5	Max. device pressure test	36
10.6	Gas Flow Rate Test	36
11	Electromagnetic Compatibility.....	37
11.1	Impact of Mobile and Portable HF Communication Devices.....	37
11.2	Electrical Connections.....	37
11.3	Guidelines and Manufacturer's Statement – Electromagnetic Emissions	38
11.4	Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity.....	39
11.5	Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity - for the Device	
	AirSeal® i.F.S.	40
11.6	Recommended Safety Distances Between Portable and Mobile RF Telecommunications Devices and the	
	AirSeal® i.F.S.	41
12	Informational, Warning and Error Messages	42
13	Technical Data	44
14	AirSeal®iFS Warranty.....	46

15	Test Log.....	48
15.1	Test Log.....	48
15.2	Return Form	49
	Index	50

1 Important User Notes

Read the manual carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing staff or service personnel, or
- to damage or malfunction of device and/or accessories.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the supplied product through continued product development.

The words WARNING, CAUTION, and NOTE carry special meanings. Sections marked with these words must be read especially attentively.

Subject to technical changes

Please note

WARNING!
The safety and/or health of the patient, user, or a third party are at risk. Comply with this warning to avoid injury to the patient, user, or third parties.

CAUTION!
These paragraphs include information provided to the operator concerning the intended and proper use of the device or accessories.

NOTE!
Here you will read information about the maintenance of the device or the accessories.

Federal Law (only for U.S. market)

Exclusion of liability

Authorized trained personnel

Intended use

Care and maintenance

Contamination

Waste management

2 Safety Information

U.S. federal law restricts use of this device to use by or on the order of a physician.

The manufacturer is not liable for direct or consequential damage and the warranty is null and void if:

- the device and/or the accessories are improperly used, prepared, or maintained,
- the instructions and rules in the manual are not adhered to,
- non-authorized persons perform repairs, adjustments, or alterations on or to the device or accessories,
- non-authorized persons open the device,
- the prescribed inspection and maintenance schedules are not adhered to.

Receipt of technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on or to the device or accessories.

Only a SurgiQuest technician may perform repairs, adjustments, or alterations on the device or accessories and use the service menu. Any violation will void the manufacturer's warranty.

The device may be used only as intended.

The service and maintenance of the device and its accessories has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is properly connected and functional.

Before shipping, decontaminate device and accessories in order to protect the service personnel. Follow the instructions listed in this manual. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to reject contaminated products for repair.

This symbol indicates that the waste of electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately instead. Please contact the manufacturer or an accordingly authorized disposal or waste management company for further information.

2.1 Hazards

WARNING!
Condensation / Water penetration
Protect device from moisture. Do not use if moisture has penetrated the device.

WARNING!
Technique and procedures
Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.

WARNING!

Check all factory settings.

Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.

**WARNING!**

Original accessories

For your own safety and that of your patient, use only original accessories.

**WARNING!**

Not explosion-proof

The device is not explosion-proof. Do not use in an area where flammable anesthetic gases are present.

**WARNING!**

Risk of electrical shock

To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.

**WARNING!**

Professional qualification

This manual does not include descriptions or instructions for surgical procedures/techniques. It is not suitable for training physicians in the use of surgical techniques. Medical peripherals and devices may be used only by physicians or medical assistants with the appropriate technical/medical qualifications working under the direction and supervision of a physician.

**WARNING!**

Functional test

The functional test must be performed prior to each surgery. Because the functional test is performed during initial start up, the unit must be power cycled (off/on) prior to each surgery.

**WARNING!**

Sterile mediums and accessories

Always work exclusively with sterile substances and mediums, sterile fluids, and sterile accessories if so indicated.

**WARNING!**

Replacement device and accessories

In case the device or any of the accessories fail during surgery, a replacement device and replacement accessories should be kept within close proximity to be able to finish the operation with the replacement components.

WARNING!
Cleaning the device
Do not sterilize the device.

WARNING!
Replacing fuse
Replace the fuse only with a fuse of the same type and rating.

WARNING!
Device-inherent dangers
Read the warnings specific to this device in chapter 3.1 "Device-inherent Dangers".

CAUTION!
To avoid the risk of electrical shock, only use this device when connected to a properly grounded power supply network.

CAUTION!
Endoscope
The device may only be connected with endoscopes designed for and featuring the technical specification permitting such a combined use. Any utilized endoscopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600.

WARNING!
Reprocessing of sterile disposable products
Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product! Do not reprocess the product.

3 Device Purpose

The SurgiQuest AirSeal® i.F.S. Intelligent Flow System (hereafter referred to as AirSeal® i.F.S.) is intended for use in diagnostic and/or therapeutic endoscopic procedures to distend the peritoneal cavity by filling it with gas, to establish and maintain a path of entry for endoscopic instruments and to evacuate surgical smoke. The cannula of the AirSeal® i.F.S. is indicated for use with or without visualization.

The device should not be used to fill an abdomen with CO₂ if a laparoscopy is contraindicated. Please consult the manual of your laparoscope for absolute and relative contraindications. The device is not suitable for hysteroscopic insufflations, i.e., it may not be used to distend the uterus.

Intended use

Contraindications

3.1 Device-inherent Dangers

WARNING!

Positioning the patient

Positioning the patient lower than the device can prevent body fluids from leaking into the tube set. Actual pressure may increase and fluid may penetrate the insufflation tube if the patient is repositioned during surgery.

WARNING!

Removing the insufflation tube

When insufflation is no longer required, press "stop" and disconnect the insufflation tubing when appropriate.

WARNING!

Backflow

Body secretions or contaminated gas may backflow into the device through the insufflation tube set if

- a filter is not used,
- the actual pressure is higher than the nominal pressure or
- the automatic venting valve is activated.

WARNING!

Gas flow

A high gas flow can occur due to large leaks within the surgical system or instrument. This can result in a false actual pressure reading, which in turn may endanger the patient. In case of a disrupted gas flow, you should therefore inspect device, tube, and instruments immediately. Surgical procedures should be performed with a gas flow of 4 to 10 l/min. An even lower gas flow is recommended for diagnostic purposes.

WARNING!

Gas supply

Maintain adequate gas supply at all times.

WARNING!

Contamination

Do not use device and/or accessories if signs of contamination are detected. Make sure the device or/and accessories can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

WARNING!
Fatigue symptoms
When there is a high level of CO₂ consumption, you should make sure to supply the operating area with enough fresh air, since an increasing CO₂ level in the air can cause the medical personnel to suffer fatigue symptoms, an inability to concentrate, unconsciousness, or even death.

WARNING!
The venting rate of the automatic venting system is limited. Always monitor the actual pressure when using additional insufflation sources.

WARNING!
Contaminated filter
Replace a contaminated filter immediately during surgery to ensure unhindered gas flow.

WARNING!
Connecting the tube
Always use the proper tube set for the device. The tube outlet may only be connected to instruments which are intended for intra-abdominal CO₂ insufflation.

WARNING!
Electronic device control
Do not close the valve at the trocar sleeve during surgery. The electronic control unit of the device adjusts the actual pressure as desired.

WARNING!
Medically pure CO₂
Make sure to use only medically pure CO₂. Other gases (i.e., helium, N₂O, argon), mixtures of gases, high pressure compressed gases, gases with entrapped liquids, or polluted gases must not be used with this device.

WARNING!
Service connection
Connected devices have to comply with the EN 60950 standard. Do not connect a device to the service connection during surgery.

CAUTION!
Electrical interference
This device was tested for electrical safety and electromechanical compatibility. However, if you detect or suspect such interference, please follow these suggestions:

- Move this, the other or both devices to a different location
- Increase distance between used devices
- Consult an electro-medical expert

WARNING!**Peripheral devices**

Additional peripheral equipment connected to interfaces of the AirSeal® i.F.S. must meet the requirements of the following specifications: IEC 60601-2-18 / EN 60601-2-18 for endoscopic devices and IEC 60601-1 / EN 60601-1 for electrical medical devices. All configurations have to comply with IEC 60601-1 / EN 60601-1 specifications. Whoever connects additional equipment to signal output or signal input is considered the system configurator and as such is responsible for complying with requirements of the standard IEC 60601-1 / EN 60601-1.

**WARNING!****Idiosyncratic reactions**

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO₂ absorption (idiosyncratic reaction).

**WARNING!****CO₂ absorption**

CO₂ is absorbed during insufflation (intravasation). This means the body absorbs part of the CO₂ gas used for insufflation. CO₂ concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO₂. High pressure or a high gas flow promotes CO₂ absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the max. intra-abdominal pressure of 30 mmHg.

**WARNING!**

The insufflation of CO₂ should be done carefully and while monitoring the patient's response. The user, particularly the anesthetist, should be informed about possible cardiovascular and respiratory problems of the patient and monitor these intra-operatively.

**WARNING!****Metabolic and cardiac reactions**

Insufflating CO₂ may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**WARNING!****Hypothermia/monitoring body temperature**

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. To reduce this risk, minimize high gas flows due to large leaks, and the use of cold irrigation and infusion solutions. Always monitor the patient's body temperature during the entire surgical procedure.

WARNING!

Lowest flow and pressure

Depending on age and health condition of the patient, the smallest possible flow and pressure for establishing the pneumoperitoneum should be selected.

WARNING!

Laparoscopy with CO₂ should not be performed on children with cardiovascular problems.

WARNING!

Dehydration

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

WARNING!

Embolism/insufflation of internal organs

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in gas embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. gas embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

WARNING!

Emphysema

Incorrect placement of a cannula or a trocar into subcutaneous tissue may lead to emphysema. To reduce the risk, use a low gas flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Long surgeries (> 200 min.), the use of many access points, duration and size of leaks at these points may also contribute to emphysema. Be sure to close leakages in trocar access points immediately.

WARNING!

Additional insufflation sources

The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

4 Initial Device Startup

Always check all parts and accessories of the device immediately after receiving the shipment. The manufacturer considers only replacement claims that have been immediately submitted or reported to a sales representative or an authorized agent.

Place the device on a flat surface and install in dry environment. The ambient temperature and humidity must meet the requirements mentioned in chapter 13 "Technical Data", page 44.

Ensure good ventilation of the device. The vents are located on the underside and on the back (see 5.2 "Device Rear" and 5.3 "Bottom of the Device").

WARNING!

Not explosion-proof

The device is not explosion-proof. Do not use in an area where flammable anesthetic gases are present.

Delivery inspection

Setting up the device

Mains connection

CAUTION!

Check to make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.

Make sure the connection data and technical specifications of the power supply comply with DIN VDE or national requirements. The mains power supply cable must be plugged into a properly installed safety wall plug (see DIN VDE 0100-710). Read the device label located in rear of device (type plate) to determine the operating voltage of the device.

The power connection must be equipped with a grounding contact. Use the original power cable (if included in scope of delivery) to establish a connection between the mains wall socket and the non-heating device plug located in the rear of the device.

Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 or IEC 60320. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.

Integrate the device into the potential equalization system as specified by local safety rules and regulations.

Grounding contact

Only for U.S. operators

Potential equalization

4.1 Gas Connection


WARNING!

Medically pure CO₂

Make sure to use only medically pure CO₂. Other gases (i.e., helium, N₂O, argon), mixtures of gases, high pressure compressed gases, gases with entrapped liquids, or polluted gases must not be used with this device.

Use a high-pressure tube to connect a CO₂ gas cylinder to the rear gas inlet connection or connect to centralized CO₂ gas supply.

4.1.1 Connecting a Gas Cylinder

CAUTION!
Always use a high-pressure tube to connect gas cylinder and device.

The gas cylinder must be in a vertical position. The gas bottle pressure may not exceed 80 bar/1160 psi. For the device start, the gas cylinder pressure must be at least 25 bar/363 psi.

CAUTION!
Gas cylinders with riser pipe can release dirt and oily fluids into the device. Do not use a gas cylinders with riser pipe.

4.1.2 Connecting to Central Gas Supply

Use the following device connectors and high-pressure tubes available as additional equipment to connect to a central gas supply (house supply):

- Conector and hose for house gas supply NIST or
 - Connector and hose for house gas supply DISS
1. Attach the high-pressure tube to the gas connection.
 2. Fasten the high-pressure tube with the nut.
 3. Tighten the nut.

NOTE!
The default gas settings are selected at the factory and should be changed as necessary by the user. If a gas supply of greater than 15 bars/218 psi is supplied to the device, regardless if the default setting is set to bottle gas or house gas, the device will automatically default to bottle gas mode.

Gas supply displays

4.1.3 Gas Consumption Display

The gas consumption display indicates the insufflated volume of CO₂ in liters since the last resetting of the display. The display depicts values between 0 and 999 liters.

The gas consumption display can be reset by pressing the **RESET** key and then returns to 0.

The status of the gas supply is monitored by the device and indicated with symbols and acoustic signals.

The following gas bottle pressures are displayed:

Gas supply with gas bottle

	> 40 bar/580 psi
	30 - 40 bar/435 - 580 psi
	15 - 30 bar/218-435 psi;
	5 - 15 bar/73-218 psi; warning screen "Gas level low. Prepare to change gas bottle" appears and audible signals can be heard. If insufflation is stopped, the warning message "Change gas bottle!" is displayed and insufflation cannot be restarted until the pressure is > 15 bar/218 psi.
	< 5 bar/73 psi; If insufflation is started, the warning message "Change gas bottle!" is displayed and acoustic signals (beeps) are emitted. The gas bottle should be changed immediately. If insufflation is stopped, the warning message "Change gas bottle!" is displayed and insufflation cannot be restarted. Smoke evacuation level will switch to LOW until tank is replaced. While in AirSeal mode, a countdown of 100 s is displayed during which the empty gas bottle can be changed while maintaining abdominal pressure.
	0 bar/0 psi; If insufflation is started, the warning message "Change gas bottle! Insufflation stopped!" is displayed and acoustic signals (beeps) are emitted. The gas bottle should be changed immediately. If insufflation is stopped, the warning message "Change gas bottle!" is displayed and insufflation cannot be restarted.

The following house gas supply pressures are displayed:

House gas supply

House gas supply pressure OK

House gas supply pressure too low


5 Operating the Device

5.1 Device Front

Familiarize yourself with the control and function elements at the front of the device.

Fig. 5-1 Front of the Device

- ① ON/OFF switch
- ② Touch screen display
- ③ Receptacle for AirSeal and Smoke Evacuation mode
- ④ Lever for locking tube set in place
- ⑤ Insufflation tube connection for Standard Insufflation Mode (barbed connection)


5.2 Device Rear

Familiarize yourself with the connection elements at the rear of the device.

Fig. 5-2 Rear of the Device


- ⑥ USB port (service only)
- ⑦ Connection for potential equalization
- ⑧ Gas connection
- ⑨ Ventilation slots (air outlet)
- ⑩ Type plate
- ⑪ Device data plate
- ⑫ Device connector plug with fuse holder

5.3 Bottom of the Device

Fig. 5-3 Bottom of the device

- ⑬ Ventilation slots (air intake)

WARNING!

The device is equipped with a powerful ventilation system with air intake located at the bottom of the device. The suction can be so strong that it can trap sheets of paper or soiling if they come in proximity. Keep ventilation slots free of obstructions or soiling to ensure optimal cooling of the device.

5.4 Display

EN


Fig. 5-4 Display

- ① Nominal gas flow rate level 1
- ② Actual flow rate display
- ③ Nominal gas flow rate level 2
- ④ Nominal gas flow rate level 3
- ⑤ Decreasing nominal pressure
- ⑥ Actual pressure display
- ⑦ Nominal pressure display
- ⑧ Increasing nominal pressure
- ⑨ Gas supply display
- ⑩ Gas Consumption Display
- ⑪ Reset key for consumption display
- ⑫ Smoke gas evacuation level LOW
- ⑬ Smoke gas evacuation status
- ⑭ Smoke gas evacuation level HIGH
- ⑮ Status display/error & warning messages
- ⑯ Mode
- ⑰ START/STOP key
- ⑱ MENU key
- ⑲ HOME key

The above depicted display shows all displays and keys. Additional explanations for individual elements are presented in the subsequent respective control element descriptions.

Not all displays and functions are available for all modes.

5.5 Switch Device On

1. Plug the device into the power outlet.
2. Connect the gas supply to the gas connection port and open the gas supply.
3. Make sure **no** tube set is connected before switching the device on.
4. Press the **ON/OFF** switch (see Fig. 5-1 "Front of the Device" ①). The device now initializes and then runs the initial self test. The start screen and a progress bar are depicted on the display during the initial self test.


If the initial self test was unsuccessful, an error message and information about how to possibly remedy the problem are displayed. An acoustic warning signal is emitted (see 12 "Informational, Warning and Error Messages", page 42).

5. The message **Device ready** is shown after the successful initial self test. An acoustic signal is emitted and the display depicts **Mode Selection**.

NOTE!

After the successful initial self test, factory-new devices will request the user select a language. Press the desired language.

6. Press the respective key to choose the desired operating mode (e.g. **Smoke**

Evacuation Mode).

Press the **Information symbol** for a short instruction about the desired operating mode. Tutorials for the three insufflation operating modes are available.

Press **HOME** key to return to **Mode Selection**.

5.5.1 Selecting Operating Mode

Press the corresponding key in **Mode Selection** to select the desired operating mode.

The device can be used with three different insufflation operating modes:





- 1. **Standard Insufflation Mode:** Traditional insufflation with a **Single Lumen Filtered Tube Set** and conventional cannula.
- 2. **Smoke Evacuation Mode:** Traditional insufflation and smoke evacuation with a **Bifurcated Smoke Evac Filtered Tube Set** for two conventional cannulas.
- 3. **AirSeal Mode:** Insufflation with **Tri-Lumen Filtered Tube Set** and AirSeal® cannula.

5.5.2 Insufflation Tube Sets

Please use the indicated accessories when working in the respective modes.

A conventional tube for **Standard Insufflation Mode** can be connected to the front of the device (see Fig. 5-1 "Front of the Device", page 14 (5)). The tube sets for the operating modes **Smoke Evacuation** and **AirSeal** are connected to the filter connection (see Fig. 5-1 "Front of the Device", page 14 (3)).

Single Lumen Filtered Tube Set (hereafter referred to as single lumen tube)	• Disposable conventional single lumen tube with standard insufflation filter	
Bifurcated Smoke Evac Filtered Tube Set (hereafter referred to as bifurcated tube)	Disposable insufflation tube set with filter • Double-lumen bifurcated tube, transparent tube for insufflation, blue tube for smoke gas evacuation	
Tri-Lumen Filtered Tube Set (hereafter referred to as tri-lumen tube)	Disposable insufflation tube set with filter for use with: • Triple-lumen tube only for use with AirSeal® cannula	
SurgiQuest Veress Needle Adapter Tube Extension (hereafter referred to as Veress Needle Tube Extension)	AirSeal®-Veress Needle adapter • Short section of tube for use with a Veress needle or traditional cannula with the Tri-Lumen Filtered Tube Set	

SurgiQuest AirSeal® Blade-less Optical Tip Obturator (hereafter referred to as Optical Obturator)		
SurgiQuest AirSeal® Blunt Tip Obturator (hereafter referred to as Blunt Tip Obturator)		
SurgiQuest AirSeal® Cannula (here after referred to as AirSeal® Cannula)		
AirSeal® Access Port and Low Profile Obturator with Blade-less Optical Tip		

After pressing the operating mode key, an informational window opens for the required tube (only for **AirSeal** and **Smoke Evacuation Mode**). Insert the tube set by positioning the filter canister into the AirSeal® i.F.S. and pulling down the lever (**AirSeal** and **Smoke Evacuation Mode**) or by pushing the open side of the tube set over the barbed connector (**Standard Insufflation Mode**).

5.5.3 Starting/Stopping Insufflation

Start insufflation:

Press **START** to start insufflation.

Activated insufflation:

The status line depicts **Initial Insufflation**. The **Initial Insufflation** phase remains in operation until the set nominal pressure has been reached through the insufflation for the first time.


CAUTION!


For the safety of the patient please fill the tube set with CO₂ gas prior to beginning the insufflation by activating the insufflation for a few seconds and then turning it off again before introducing the insufflation instrument to the abdomen and beginning the surgery.

Stop insufflation:

Press **STOP** to stop insufflation.

EN

5.5.4 Power Device Off

Use the **ON/OFF** button to turn the device off.

6 Using and Controlling the AirSeal® i.F.S. in its Different Modes

The device can be used with three different insufflation operating modes (see 5.5.2 "Insufflation Tube Sets"):

1. **Standard Insufflation Mode:** Traditional insufflation with a **Single Lumen Filtered Tube Set** and conventional cannula.
2. **Smoke Evacuation Mode:** Traditional insufflation and smoke evacuation with a **Bifurcated Smoke Evac Filtered Tube Set** for two conventional cannulas.
3. **AirSeal Mode:** Insufflation with **Tri-Lumen Filtered Tube Set** and AirSeal® cannula.

6.1 Standard Insufflation Mode

The **Standard Insufflation Mode** is intended for use with a **Single Lumen Filtered Tube Set** and standard cannulas.

1. Switch device on (consult section 5.5 "Switch Device On", page 15 for additional information).
2. Press the **Standard Insufflation Mode** key to select this operating mode.
3. Connect the Single Lumen Filtered Tube Set.

Inserting the Single Lumen Filtered Tube Set

To be carried out by non-sterile technician:



- Open packaging of the Single Lumen Filtered Tube Set.
- Have a sterile technician remove the tube set from the inside of the package.

To be carried out by sterile technician:

- Keep the distal tube Luer lock connector in the sterile area and hand the proximal tube end with filter to the non-sterile technician.
- Connect the Luer lock connector with the instrument (e.g. inflow cannula). Open inflow valve at instrument.

To be carried out by non-sterile technician:

- Attach the insufflation tube set to the insufflation tube set connection at the front of the device (barbed connector).

4. Choose desired settings.

Setting the nominal flow and nominal pressure (increase/decrease) is possible during insufflation or while insufflation is stopped.

Setting nominal flow:

The device features three different flow level rates that can be adjusted in the user menu (see chapter 8 "User Menu" for additional information):

Level 1 -> 5 l/min

Level 2 -> 20 l/min

Level 3 -> 40 l/min

- To set the flow rate, press the key (1), (3) or (4).

Setting the nominal pressure:

Press the key ▼ or ▲ ((5) or (8)) to set the nominal pressure. Values may range from **5 to max. 20 mmHg** and can be adjusted in increments of 1. The

starting value for the pressure can be set in the user menu within a range of 5 - 15 mmHg.

- Scrolling is enabled by keeping the ▼ or ▲ key depressed longer than 1.5 seconds.

Safety limit:

When wanting to **Increase** the nominal pressure to > 15 mmHg, the status line depicts the message **Safety Limit: > 15 mmHg**. This is where the recommended range for the intra-abdominal pressure ends. Pressing the nominal pressure key ▲ again does not increase the pressure any further. Release the key for 2 seconds. Then a value up to 20 mmHg can be set.

CAUTION!

Exceeding this safety limit is to be decided by and the responsibility of the user/operator.

5. Start insufflation with the **START** key.

The status line depicts **Initial Insufflation**. The initial insufflation ends when the nominal pressure value has been reached for the first time. After the initial insufflation, the flow rate automatically increases to the nominal gas flow specified for level 3. The status line depicts **Standard Insufflation Active**.

NOTE!

For the safety of the patient, it is advised to start insufflation with a Veress needle and at the lowest flow rate (level 1, key ①). Exchange the Veress needle with a conventional cannula once the nominal pressure is reached and the initial insufflation phase ends.

CAUTION!

Always use Single Lumen Filtered Tube Set with filter (see 5.5.2 "Insufflation Tube Sets").

WARNING!

Positioning the patient

Always position the patient lower than the device to prevent body fluids from leaking into the insufflation tube. Actual pressure may increase and fluid may penetrate the insufflation tube if the patient is repositioned during surgery. If this occurs, immediately disconnect the insufflation tube. When the patient is repositioned onto his or her side, internal tissue may block the insufflation channel. Always insufflate through the elevated side of the patient.

WARNING!

Backflow

Body secretions or contaminated gas may backflow into the device through the insufflation tube if

- a filter is not used,
- the actual pressure is higher than the nominal pressure or
- the automatic venting valve is activated.

6. Stop insufflation with the **STOP** key.

6.2 Smoke Evacuation Mode

The **Smoke Evacuation Mode** is intended for standard insufflation with continuous evacuation of smoke from the abdomen as it may occur during the use of standard ultrasonic, laser and electric energy devices. Use the Bifurcated Smoke Evac Filtered Tube Set with standard cannulas when working in this mode. The **Smoke Evacuation Mode** is automatically activated when the second line of the bifurcated tube set (blue-colored tube) is attached to an open second cannula.

1. Switch device on.
2. Press the **Smoke Evacuation Mode** key.
3. Insert the Bifurcated Smoke Evac Filtered Tube Set.

Inserting the Bifurcated Smoke Evac Filtered Tube Set

To be carried out by non-sterile technician:


- Open packaging of the Bifurcated Smoke Evac Filtered Tube Set.
- A sterile technician then removes the tube set from the inside of the package.

To be carried out by sterile technician:

- Keep the Luer lock connectors of the distal bifurcated tube in the sterile area and hand the proximal tube end connected to the filter to the non-sterile technician.
- Connect the clear Luer lock connector of the bifurcated tube with the conventional cannula. Open the inflow valve.
- Connect the blue Luer lock connector of the bifurcated tube to a different conventional cannula to activate **smoke evacuation**.

To be carried out by non-sterile technician:

- Insert the Smoke Evacuation filter housing into the receptacle at the front of the device. Use the lever to lock the filter housing in place or to unlock and release it.

4. Choose desired settings.

Setting the nominal flow and nominal pressure (increase/decrease) is possible during insufflation or while insufflation is stopped.

Setting nominal flow:

The device features three different flow level rates that can be adjusted in the user menu (see chapter 8 "User Menu" for additional information):

Level 1 -> 5 l/min

Level 2 -> 20 l/min


Level 3 -> 40 l/min

- To set the flow rate, press the key (1), (3) or (4).

Setting the nominal pressure:

Press the key ▼ or ▲ ((5) or (8)) to set the nominal pressure. Values may range from **5 to max. 20 mmHg** and can be adjusted in increments of 1. The starting value for the pressure can be set in the user menu within a range of 5 - 15 mmHg.

- Scrolling is enabled by keeping the ▼ or ▲ key depressed longer than 1.5 seconds.

Safety limit:

When wanting to **Increase** the nominal pressure to > 15 mmHg, the status line depicts the message **Safety Limit:> 15 mmHg**. This is where the recommended range for the intra-abdominal pressure ends. Pressing the nominal pressure key ▲ again does not increase the pressure any further. Release the key for 2 seconds. Then a value up to 20 mmHg can be set.

CAUTION!

Exceeding this safety limit is to be decided by and the responsibility of the user/operator.

5. Start insufflation by pressing the **START** key.

The status line depicts **Initial Insufflation**. The initial insufflation ends when the nominal pressure value has been reached for the first time. After the initial insufflation, the flow rate automatically increases to the nominal gas flow specified for level 3. The device prompts to connect the second line to enable the Smoke Evacuation function.

To enable smoke gas evacuation, connect the blue tube of the bifurcated tube to a second cannula.

NOTE!

For the safety of the patient, it is advised to start insufflation with a Veress needle and at the lowest flow rate (level 1, key ①). Exchange the Veress needle with a conventional cannula once the nominal pressure is reached and the initial insufflation phase ends.

6. Stop insufflation by pressing the **STOP** key.

When working in the **Smoke Evacuation Mode**, it is possible that body fluids are extracted from the surgical field over the evacuation line. In order to prevent contamination of the device, these are retained by the fluid trap of the Smoke Evacuation tube set filter component. The capacity of the fluid trap is limited. A warning message and an audible signal will be emitted if the fill level **LOW** is reached. Cannula and tubing placement should be checked to make sure no more fluid enters the filter. At this point, change the filtered tube set. Insufflation can be continued with full functionality. If fluid level **HIGH** is reached, the Smoke Evacuation function will shut down but insufflation is still active.

6.3 AirSeal Mode

The **AirSeal Mode** is intended to establish and maintain a pathway of entry for laparoendoscopic instruments to the field of surgery within laparoscopic procedures. The abdominal pressure is maintained by an air barrier within the cannula system that allows valve-free access to the body cavity and that does not collapse when introducing instruments.

WARNING!

Please read the respective instruction manual before using the AirSeal® cannula.

6.3.1 AirSeal® Optical Obturator and Cannula System

The AirSeal® Optical Obturator and Cannula System is packaged in three configurations:

- Optical obturator


Fig. 6-1 AirSeal® Optical Obturator

- ① Bladeless optical tip obturator
- ② Cannula
- ③ Manifold plug
- ④ Cannula latch

- Blunt obturator with suture anchor


Fig. 6-2 AirSeal® Blunt Obturator with Suture Tie Down

- ① Blunt obturator
- ② Cannula
- ③ Manifold plug
- ④ Spring anchor

- Access Port and Low Profile Obturator with Bladeless Optical Tip


Fig. 6-3 AirSeal® Access Port and Low Profile Obturator with Bladeless Optical Tip

- ① Bladeless optical tip obturator
- ② Cannula
- ③ Manifold plug
- ④ Cannula latch

The AirSeal® Optical Obturator and Cannula System is composed of a sterile single patient use instrument consisting of an obturator ① and a cannula ② (see Fig. 6-1, Fig. 6-2, and Fig. 6-3). The obturator may be used with or without visualization for primary and secondary insertions.

If the Optical cannula is being used for primary entry, use of a laparoscope to provide visualization to enhance safe abdominal entry is advised.

1. Use sterile techniques to remove the instrument from the package.
2. The obturator ① and the cannula ② (see Fig. 6-1, Fig. 6-2, and Fig. 6-3) are packaged together.
3. **When using optical obturator**
 - Connect an appropriate 0° endoscope to the light supply and monitor (see manufacturer's instructions).
 - Insert the endoscope into the opening at the proximal end of the obturator until it reaches the distal tip of the obturator. Secure the endoscope in the obturator using the scope lock.

6.3.2 Initial Insufflation

When working in the **AirSeal Mode**, it is possible to initiate insufflation either through the open cannula, through the cannula with inserted obturator or through a Veress needle.

CAUTION!

The manufacturer strongly advises to begin the insufflation by using the Veress needle

Initial insufflation with Veress needle entry

1. Switch device on (consult section 5.5 "Switch Device On", page 15 for additional information).
2. At the conclusion of the initial self test, press the **AirSeal Mode** key to select this operating mode.
3. Insert the Tri-Lumen Filtered Tube Set.

Inserting the Tri-Lumen Filtered Tube Set

To be carried out by non-sterile technician:

- Open packaging of the Tri-Lumen Filtered Tube Set.
- Have a sterile technician remove the tube set from the inside of the package.

To be carried out by sterile technician:

- Keep the distal end of the Tri-Lumen Filtered Tube Set AirSeal® cannula connector in the sterile area and hand the tube end with filter to the non-sterile technician.

Option a) Initial insufflation with Veress needle or conventional cannula (advised)


Use the Veress Needle adapter to connect the tri-lumen tube with the Luer lock connector of a conventional Veress needle or a conventional cannula (see 5.5.2 "Insufflation Tube Sets"). The Veress adapter comes connected to the tri-lumen tube upon delivery.

Option b) Initial insufflation with AirSeal® cannula

Connect the Tri-Lumen Filtered Tube Set with the AirSeal® cannula (see 6.3.1 "AirSeal® Optical Obturator and Cannula System").

To be carried out by non-sterile technician:

- Insert the Tri-Lumen Filtered Tube Set into the receptacle at the front of the device. Use the lever to lock the filter housing in place or to unlock and release it.

4. Insert the Veress needle or AirSeal® cannula in accordance with proper laparoscopic technique.
5. Choose desired settings.
Setting the nominal flow and nominal pressure (increase/decrease) is possible during insufflation or while insufflation is stopped.

Setting nominal flow:

The device features three different flow level rates that can be adjusted in the

user menu:

Level 1 -> 5 l/min

Level 2 -> 20 l/min

Level 3 -> 40 l/min

- To set the flow rate, press the key (1), (3) or (4).

Setting the nominal pressure:

Press the key ▼ or ▲ (5) or (8), to set the nominal pressure. Values may range from **5 to max. 20 mmHg** and can be adjusted in increments of 1. The starting value for the pressure can be set in the user menu within a range of 5 - 15 mmHg.

- Scrolling is enabled by keeping the ▼ or ▲ key depressed longer than 1.5 seconds.

Safety limit:

Increase the nominal pressure to > 15 mmHg, status line depicts the message **Safety Limit: > 15 mmHg**. This is where the recommended range for the intra-abdominal pressure ends. Pressing the nominal pressure key ▲ again does not increase the pressure any further. Release the key for 2 seconds. Now you can set a value up to 20 mmHg.

CAUTION!

Exceeding this safety limit is to be decided by and the responsibility of the user/operator.

6. Start insufflation with the **START** key.
7. The status line depicts **AirSeal initiating**. Initial Insufflation ends when the AirSeal® System is active.
Note that the use of a Veress needle is automatically detected by the device. The AirSeal® recirculation gas flow is **NOT ACTIVATED** when performing initial insufflation with a Veress needle in the **AirSeal Mode**.

NOTE!


For the safety of the patient, it is advised to start insufflation with a Veress needle and at the lowest flow rate (level 1, key (1)).

8.

Initial Insufflation with a Veress Needle

I. Once the set pressure is achieved (information screen on display), detach the tri-lumen tube from the Veress needle tube extension. Remove the Veress needle and insert the AirSeal® cannula with inserted obturator. Remove the manifold plug from the side of the cannula housing and attach the tri-lumen tube to the AirSeal® cannula manifold and tighten.

Push the **AirSeal** button on the device to activate the recirculation gas flow for building up the air seal. Follow on screen messages until prompted to remove obturator.

CAUTION!

Should the AirSeal Mode button not be pushed before removing the obturator from the cannula, the pressure in the abdomen can not be maintained. Make sure the recirculation gas flow is activated before removing the obturator.

Initial Insufflation with AirSeal® Cannula

II. Once the recirculation gas flow is activated, the obturator can be removed. Follow on screen messages.

9. When procedure is completed stop insufflation by pressing the **STOP** key.
10. When working in the **AirSeal Mode**, it is possible that body fluids can be extracted from the surgical field over the evacuation line and into the filter housing. In order to prevent contamination of the device, bodily fluid is retained by the fluid trap of the housing component. The capacity of the fluid trap is limited. A warning message and an audible signal will be emitted if the fill level **LOW** is reached. Cannula and tubing placement should be checked to make sure no more fluid enters the filter. At this point, change the filtered tube set. Insufflation can be continued with full functionality. If fluid level **HIGH** is reached, the AirSeal function will shut down. Insert obturator in AirSeal Cannula to maintain pressure with normal insufflation.

7 Safety Functions

7.1 General Safety Functions

The device is equipped with an automatic venting system.

When the insufflator detects that the nominal pressure was exceeded by more than 3 mmHg for longer than 3 seconds, it automatically activates the venting system. The status line then depicts **Venting active!**

WARNING!

The venting rate of the automatic venting system is limited. Always monitor the actual pressure when using additional insufflation sources.

If the overpressure of more than 3 mmHg exists longer than 5 seconds, the status line depicts a corresponding message and a warning signal is emitted.

CAUTION!

In case of overpressure above 30 mmHg for longer than 5 seconds insufflation is deactivated.

When tube, Veress needle, or AirSeal® cannula are blocked, the message **Occlusion** is depicted and a warning signal is emitted. The actual pressure display depicts **0**.

The warning signal can be deactivated in the user menu (see chapter 8 "User Menu"). Remove the occlusion if this is not done intentionally, e.g. with a closed valve. Check tubes or cannula for proper fit (not in tissue) and check for possible blockages.

7.2 Contamination Alarm

Contamination occurs when fluid has penetrated the device via the conventional barbed fitting tube connection. A contamination message and an audible signal occurs. It is possible to conclude the currently ongoing surgery with this device. Insufflation is no longer possible after turning the device off and back on using the **ON/OFF** switch. This is to prevent cross-contamination.

The display depicts **Contamination** if you are switching on an already contaminated device. The device can no longer be used.

The contaminated device has to be clearly marked as contaminated and sealed in two separate protective layers of safety foil.

Make sure the device can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.

Venting system

Overpressure

Occlusion

Contamination

Stop insufflation / System calibration

Severed tube connection

Leakage warning

Fluid level LOW (AirSeal Mode and Smoke Evacuation Mode)

Fluid level HIGH (Smoke Evacuation Mode)

Fluid level HIGH (AirSeal Mode)

Fluid in filter before surgery

7.3 Safety Functions in AirSeal Mode

The pressure is reduced slowly before the system comes to a complete stop when stopping the AirSeal® insufflation.

The stopping procedure consists of an automatic calibration routine so that the unit keeps running a few seconds after stop key is pressed. During this process no pressure or insufflation is built up.

If the connection of the AirSeal® filter tube set with the AirSeal® cannula is interrupted or severed, the system stops and the warning **Check connection of Tri-Lumen Filtered Tube Set!** is depicted.

If a large leak makes it impossible to maintain the AirSeal® function, the warning message **Check for excessive leakage!** is depicted.

7.4 Fill Level Display

Tri-lumen and the bifurcated filtered tube sets are equipped with fluid sensors to monitor the fluid level in the filters and to warn users of a possible device contamination.

When working in the **AirSeal Mode** and **Smoke Evacuation Mode**, body fluids may be flushed from the surgical field through the evacuation line and into the filter housing. In order to prevent contamination of the device, bodily fluid is retained by the fluid trap of the housing component. The capacity of the fluid trap is limited.

A warning message and an audible signal will be emitted if the fill level **LOW** is reached. Cannula and tubing placement should be checked to make sure no more fluid enters the filter. Change the filtered tube set, and insufflation can be continued with full functionality.

If fluid level **HIGH** is reached in Smoke Evacuation Mode: Smoke Evacuation function will stop but insufflation can still be continued.

If fluid level **HIGH** is reached in AirSeal Mode: AirSeal function will stop, but surgery may continue using standard insufflation.

This requires that the obturator is reinserted into the AirSeal® cannula or the tri-lumen tube is connected with the tube adapter to a conventional cannula after conventional cannula is inserted. If AirSeal® cannula is removed, the surgeon may wish to insert conventional cannulas in same entry point.

Insert a new Tri-Lumen tube to continue with the AirSeal Insufflation.

A warning screen will immediately appear if the fluid fill level in the filter housing of the tri-lumen tube sets and bifurcated tube sets is at low or high upon inserting one of those tube sets. Insufflation will not start until a new tube set is inserted.

8 User Menu

While insufflation is stopped, press the **MENU** key (Fig. A) to open the user menu (Fig B).


Fig. A

Fig. B

Device parameters can be changed in the user menu. The following pages provide an overview as well as a detail description.

Access to the service menu is restricted to trained and authorized service personnel.


Access to the service menu

8.1 Setting First Nominal Pressure

In the user menu, tap the **Starting Pressure** key to access the setting.

Select the desired **Starting Pressure** for the currently selected insufflation mode.

Press the ▲ or ▼ keys to increase or decrease the nominal pressure.


Insufflation operating mode	Factory setting	Range
AirSeal Mode	15 mmHg	5-15 mmHg
Smoke Evacuation Mode	15 mmHg	5-15 mmHg
Standard Insufflation Mode	15 mmHg	5-15 mmHg

8.2 Gas Flow Rates

In the user menu, tap the **Gas Flow Rates** key to access the selection. You can select one of three gas flow rates for each insufflation mode.

Press the ▲ or ▼ keys to increase or decrease the gas flow.

Insufflation operating mode	Range
AirSeal Mode (Initial Insufflation)	Gas flow between 1 and 40 l/min
Smoke Evacuation Mode	Gas flow between 1 and 40 l/min
Standard Insufflation Mode	Gas flow between 1 and 40 l/min

8.3 Setting Smoke Evacuation Level

In the user menu, tap the **Smoke Evacuation Level** key to access the selection. For AirSeal Mode a pre-setting can be chosen.

AirSeal Mode	LOW
	HIGH

8.4 Setting General Volume

In the user menu, tap the **Audio/Display** key to access the signal volume selection.

You can select one of four volume levels for the warning signals as well as the info and alert signals, all in the form of acoustic signals. The selection/setting applies to all insufflation operating modes.

Press the ▲ or ▼ key to set the desired volume.

Warning signals	Alert signals	Info signals
I (Low)	I (Low)	I (Low)
II (Medium)	II (Medium)	II (Medium)
III (High)	III (High)	III (High)
IIII (Very High)	IIII (Very High)	IIII (Very High)

8.5 Setting Brightness

In the user menu, tap the **Audio/Display** key to access the brightness selection.

Use the **Brightness** option to adjust the screen display to the ambient light conditions within the operating room.

8.6 Setting Audible Occlusion Signal

In the user menu, tap the **Audio/Display** key to access the audible occlusion signal selection.

The audible signal which comes with the occlusion alert can be switched off. The alert message remains nevertheless.

8.7 Setting Gas Supply

In the user menu, tap the **Gas Supply** key to access the gas supply selection.

You can select house or bottle gas supply.

House gas button will be greyed out if pressure is > 15 bar/218 psi. To select house gas supply, close and disconnect gas bottle from device.

Press the corresponding key to select the desired gas supply.

In this menu item, press the ▲ or ▼ key to select the desired volume for **No Gas Alarm**.

8.8 Setting the Language


In the user menu, tap the **Language** key to access the language selection.

Select one of ten available languages.

Language

English (factory setting)	Japanese
French	Korean
German	Portuguese
Spanish	Russian
Italian	Chinese


Press the corresponding language key to select the desired language or press the arrow key to access to further languages.

8.9 Checking Software Version


In the user menu, tap the **Software Version** key.

The software version is depicted.

8.10 Resetting or Restoring to Factory Setting

In the user menu, tap the **Reset** key. This allows you to reset all setting done in the user menu to factory default.

Authorized trained personnel

Manufacturer's specifications

Two-year maintenance interval

Authorized trained personnel

Unauthorized personnel

Liability

Technical documents

Certification

9 Care and Maintenance

Special care is necessary when servicing, maintaining, and storing the device and its accessories to maintain the functionality of the device and its accessories.

Only a SurgiQuest technician may perform repairs, adjustments, or alterations on the device or accessories and use the service menu. Any violation will void the manufacturer's warranty.

9.1 Cleaning the Device

- 1. Use the **ON/OFF** key to turn the device off.
- 2. Remove the power cable.
- 3. Wipe the surface of the device with a soft cloth moistened with an appropriate disinfectant (for example Meliseptol® rapid). The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the device.

NOTE!
Do not sterilize the device.

9.2 Annual Inspection

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. This inspection has to be carried out once a year. The tests are described in chapter 10 "Annual Inspection", page 34.

Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

9.3 Maintenance by Authorized Service Technician

An authorized service technician has to inspect and service the device at appropriate intervals to ensure the safety and functionality of the unit. The minimum service interval is two years, depending on frequency and duration of use. If the service interval is not maintained, the manufacturer does not assume any liability for the functional safety of the device.

A sticker located on the rear panel of the device will remind you of the latest date for the next service or maintenance check.

Authorized service technicians are only trained and certified by the manufacturer.

All of the service tasks, such as changes, modifications, repairs, calibrations, etc. may be carried out only by the manufacturer or manufacturer-approved trained and skilled technicians.

The manufacturer is not liable for the operational safety of the device if unauthorized persons conduct this maintenance or any other service tasks.

Unauthorized opening of the device and repairs performed by unauthorized personnel or third parties and/or changes or modifications release the manufacturer of any liability concerning the operational safety of the device.

Receiving technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on the device or accessories/peripherals.

Ask the service technician for a certificate after he or she has inspected the unit or performed any service tasks. This certificate lists the type and scope of the service as well as the date and name of the servicing company together with the signature of the service technician.

9.4 Replacing the Fuse

CAUTION!

Before replacing the fuse, check the values of the fuse to be inserted acc. to chapter 13 "Technical Data", page 44.

The fuse may be defective and is in need of replacement if:

- displays and LEDs (if available on your equipment) do not light up,
- the device does not function.

Check to make sure

- the main power supply cable is properly connected to the power supply input and to a safety socket,
- the house power supply fuse is functioning.

WARNING!

Unplug the power cable from the device before checking the fuse.

The device does **not** have to be opened to replace the fuse.

1. Switch device off.
2. Disconnect device from power supply.
3. Remove power connection cable from mains socket.
4. The fuse holder is located next to the mains socket.
5. Remove fuse holder as depicted in Fig. 9-1 "Opening the fuse holder".
6. **A** Undo the latch of the fuse holder with a small screwdriver.
7. **B** Remove the fuse holder.
8. **C** Check fuse.
9. Insert a new fuse. Use only the specified type of fuse (see chapter 13 "Technical Data", page 44).
10. Insert the fuse holder until it can be heard snapping into place.
11. Use the power cable to reconnect the shockproof safety socket with the rear mains socket.


Fig. 9-1 Opening the fuse holder

Measured values and tolerances

10 Annual Inspection

Each test conducted has to be documented with date and signature on the test log.

The following measuring tools and resources were used by the manufacturer to determine the listed measurements and tolerances:

Manometer	Range 0-100 mmHg, error class $\pm 1.6\%$
Syringe	60 ml
Silicone tube	8 mm x 2 m
T adapter	8-8-8 mm
Veress needle	length 100 mm opening diameter 1.4 mm, inner cannula diameter 1.6 mm

An authorized service technician must check the device if the specified parameters and tolerances are exceeded.

10.1 Electrical Safety Test

1. Perform a visual inspection. Make sure that
 - the fuse corresponds with the specifications indicated by the manufacturer,
 - labels and stickers on device are legible,
 - the mechanical condition of the device allows for its safe use,
 - the device is clean to ensure proper and safe functionality.
2. Perform the measurement of the ground leakage current (max. 500 μA) and contact current (max. 100 μA in normal state and max. 500 μA on first error) according to IEC 60601-1 / EN 60601-1.
3. Measure protective conductor resistance according to IEC 60601-1 / EN 60601-1. The protective conductor resistance is measured while the device is connected to the power supply. The max. value is 0.2 Ω .

As an alternative, perform safety test according to DIN EN 62353.


10.2 Basic Function Test

1. Conduct this test without the Single Lumen Filtered Tube Set connected to the AirSeal® i.F.S.
2. Use the **ON/OFF** switch to turn the device on. The device now conducts a self-test. A short acoustic signal can be heard. Press the **Standard Insufflation Mode** key to select this operating mode.
3. The factory default settings are 15 mmHg for the nominal pressure and 5 l/min for the nominal flow.
4. The following values are displayed:
 - Nominal pressure 15 mmHg**
 - Nominal gas flow rate 5 l/min [button 1]**
 - Actual pressure 0 mmHg**
 - Gas consumption 0 l**
5. Start insufflation:
 - Press the **START** key.
 - The following values are displayed: Actual pressure 0 mmHg
 - Initial insufflation** is depicted.
 - Streaming gas can be heard at the insufflation tube connection.

6. Select the max. nominal gas flow.
The following values are displayed:
Actual gas flow equals nominal gas flow rate
Actual pressure 0 mmHg
Initial Insufflation is displayed.
Streaming gas can be heard at the insufflation tube connection (barbed connector).
7. Stop insufflation:
Press the **STOP** key.
The following values are displayed:
Actual pressure 0 mmHg
Gas consumption > 0.0 l
8. Press the **RESET** key.
Gas consumption 0.0 l

The basic function check of the device is complete.

10.3 Pressure Sensor Test

1. Select the **Standard Insufflation Mode** as the operating mode.
2. Do not press the **START/STOP** key.

CAUTION!

Never use the syringe to extract gas from the device.


3. Connect a manometer (1) and an air-filled syringe (3) to the insufflation tube connection (2).
4. Use the syringe to generate a pressure of at least 10 mmHg, which is indicated on the manometer.
Actual pressure display: 10 ± 2 [mmHg]
5. Use the syringe to generate a pressure of at least 20 mmHg, which is indicated on the manometer.
Actual pressure display: 20 ± 2 [mmHg]
6. Use the syringe to generate a pressure of at least 30 mmHg, which registers on the manometer.
Actual pressure display: 30 ± 2 [mmHg]

10.4 Pressure Monitoring Test

1. Select the **Standard Insufflation Mode** as the operating mode.
2. Select a nominal pressure of 15 mmHg and adjust in the user menu the gas flow rate 1.
3. Use the syringe to generate a pressure of at least 19 mmHg, which registers on the manometer.
An acoustic warning sound is emitted with a pressure of more than 19 mmHg (for 5 seconds) and the display depicts **Overpressure**.
4. Reduce the pressure.
The warning ends when the pressure falls below 19 mmHg (nominal pressure plus 4 mmHg).

10.5 Max. device pressure test

1. Select the **Standard Insufflation Mode** as the operating mode.
2. Select the max. nominal gas flow.
3. Connect a manometer (1) and an open Veress cannula (3) to the insufflation tube connection (2).
4. Select the max. gas flow.
5. Start insufflation:
Press the **START** key.
The manometer registers a pulsing pressure increase. When the pressure stabilizes, the manometer indicates a maximum pressure between 65 and 75 mmHg.
6. Stop insufflation:
Press the **STOP** key.

10.6 Gas Flow Rate Test

Test setup with open connection, without connected insufflation tube.

1. Adjust in the user menu the gas flow rate 2 (20 liters).
2. Start insufflation: Press the **START** key.
3. Press the **RESET** key (0.0 l must be depicted).
Now start the measuring for one minute.
4. Stop the insufflation after one minute. Press the **STOP** key.
The gas consumption should be at least 17.5 - 22.5 liters.

11 Electromagnetic Compatibility

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

Precautionary measures

This device is to be used only for the purposes described in the manual and has to be installed, set up, and operated in compliance with the EMC notes and instructions.

11.1 Impact of Mobile and Portable HF Communication Devices

The emission of high frequency energy by mobile communication devices may affect the function of the electrical medical device. Operating such devices (e.g., cell phones, GPS phones) in the proximity of the electrical medical device is not recommended.

11.2 Electrical Connections

Do not touch electrical connections identified with this warning label. Do not establish a connection between these plugs and sockets without first implementing precautionary ESD (electrostatic discharge) measures.

The following are ESD precautionary measures:

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

ESD (Electrostatic Discharge) precautionary measures

Hospital employees should be informed about and trained in ESD precautionary measures.

11.3 Guidelines and Manufacturer’s Statement – Electromagnetic Emissions

The device AirSeal® i.F.S. is intended for use in the electromagnetic environment specified below. The user/operator of the AirSeal® i.F.S. should make sure the device is operated within such an environment.

Emitted interference measurements	Compliance	Electromagnetic environment guidelines
HF emission according to CISPR 11	Group 1	The device AirSeal® i.F.S. uses HF energy solely for its internal functions. Therefore, the camera's HF emission is very low and it is unlikely that devices in close proximity will experience interference.
HF emission according to CISPR 11	Class B	The device AirSeal® i.F.S. is suitable for use in all facilities including those in residential areas and those directly connected to a public utility network supplying buildings used for residential purposes as well.
Emission of harmonic oscillations according to IEC 61000-3-2	Class A	
Emission of voltage fluctuations / flickers according to IEC 61000-3-3	In compliance	

11.4 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity

The device AirSeal® i.F.S. is intended for use in an electromagnetic environment as described below. The user/operator of the AirSeal® i.F.S. should make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	Test level	Compliance	Electromagnetic environment guidelines
Discharge of static electricity (ESD) according to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	In compliance	Floors should be made from wood or concrete or covered with ceramic tiles. If the floor covering consists of synthetic material, the relative humidity should be at least 30%.
Electrical fast transients / bursts according to IEC 61000-4-4	± 2 kV for AC power lines ± 1 kV for input and output lines	In compliance	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Surges according to IEC 61000-4-5	± 1 kV normal mode voltage, ± 2 kV common mode voltage	In compliance	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Blackouts, brown-outs, and fluctuations of the power supply according to IEC 61000-4-11	< 5% U_T^* (> 95% dip in the U_T) for ½ cycle	In compliance	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment. If the user/operator of device requires the continuation of functionality after power interruptions/disruptions, it is recommended to supply the device with power from an uninterruptible power supply.
	40% U_T (60% dip in the U_T) for 5 cycles.		
	70% U_T (30% dip in the U_T) for 25 cycles.		
	< 5% U_T (> 95% dip in the U_T) for 5 s		
Supply frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m	In compliance	Magnetic fields of the mains power frequency should comply with the typical values of business and hospital environments.

*Note: U_T is the mains alternating voltage before applying the test levels.

11.5 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity - for the Device AirSeal® i.F.S.

Electromagnetic interference immunity tests	Test level	Compliance	Electromagnetic environment guidelines
Conducted HF interference quantities according to IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	In compliance	Portable and mobile wireless devices should not be used in closer proximity to the device AirSeal® i.F.S. (including cables/lines) than the recommended safety distance calculated based on the transmitting frequency and the applicable formula. Recommended safety distance: d = 1.2√P for 150 KHz to 80 MHz d = 1.2√P for 80 MHz to 800 MHz d = 2.3√P for 800 MHz to 2.5 GHz With P as the rated output of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d as recommended safety distance in meters [m]. The field strength of stationary transmitters for all frequencies tested on site ^a should be lower than the concordance level. ^b Interference is possible in the proximity of devices featuring the following pictograph.
Radiated HF interference quantities according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	In compliance	

Note 1: The higher frequency range applies for 80 and 800 MHz.

Note 2: These guidelines are probably not realizable in all cases. The distribution and spread of electromagnetic quantities differs depending on the absorption and reflection of buildings, objects, and people.

^a The field strength of stationary transmitters such as base stations of wireless phones and cell phones, ham radio operators, AM and FM radio and TV stations can theoretically not always determined in advance. A study of the installation site should be considered to determine the electromagnetic environment concerning the stationary transmitter. If the measured field strength at the proposed AirSeal® i.F.S. installation and AirSeal® i.F.S. operation site exceeds the concordance levels listed above, the AirSeal® i.F.S. should be monitored to document proper functionality and operation as intended. If unusual performance characteristics are observed, additional measures may be required such as changing orientation or the location of the device AirSeal® i.F.S..

^b The field strength should be less than 3 V/m for the frequency range of 150 kHz to 80 MHz.

11.6 Recommended Safety Distances Between Portable and Mobile RF Telecommunications Devices and the AirSeal® i.F.S.

Recommended safety distances between portable and mobile RF telecommunications devices and the AirSeal® i.F.S.			
The AirSeal® i.F.S. is intended for use in an electromagnetic environment where HF interferences are controlled. The user/operator of the AirSeal® i.F.S. can contribute to lowering electromagnetic emissions by complying with the minimum distance between portable and mobile HF telecommunications devices (transmitters) and the AirSeal® i.F.S. - depending on the output power of the communication device listed below.			
Rated output of the transmitter [W]	Safety distance based on the transmitting frequency [m]		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The safety distance d in meters [m] for transmitters with a max. rated output not listed in the table above can be calculated by applying the corresponding formula in the respective column. P is the max. rated output of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter.

Note 1: The higher frequency range applies to 80 and 800 MHz.

Note 2: These guidelines are probably not realizable in all cases. The distribution and spread of electromagnetic quantities differs depending on the absorption and reflection of buildings, objects, and people.

12 Informational, Warning and Error Messages

Status line display

Information for users displayed in the status display are depicted only temporarily and disappear again after a few seconds.

Informational, warning and error messages

Depending on the type of message, windows with green, yellow, and red frames are depicted.

Informational message	Green frame
Warning message	Orange frame
Error message	Red frame
Critical message	Red frame

Informational message	Cause
AirSeal is shutting down. System Auto Calibration in progress.	Pressure is slowly reduced. Pump keeps running for a few seconds. No pressure build up. Do not turn device off. Wait for calibration.
Connect blue tube. Smoke evacuation is initiating...	To enable the smoke gas evacuation function, connect the blue tube to a trocar and open the valve.
Remove obturator.	AirSeal function initiated.

Warning message	Cause	Remedy
Check for excessive leakage!	Long suction or leaks.	Check for possible leakage! Use suction carefully.
Fluid in Filter! Check cannula position! If fluid level increases Smoke Evacuation will stop!	Low fluid level of the fluid trap in the filter housing during Smoke Evacuation mode.	Please insert a new tube set.
Fluid in Filter! Check cannula position! If fluid level increases AirSeal insufflation will shut down.	Low fluid level of the fluid trap in the filter housing during AirSeal mode.	Please insert a new tube set.
Smoke Evacuation stopped. Insufflation active.	Leaks or low abdominal pressure.	Check for possible leakage! Use suction carefully
Smoke Evacuation stopped. Suction line blocked!	Evacuation blocked.	Check the evacuation/suction system for proper function.
Safety threshold: Pressure >15 mmHg!	The set value of 15 is reached and the user tries to increase the pressure further.	Release the key for 2 seconds. Now you can set a value up to 20 mmHg.
Starting... Do not install tube set until self test is completed!	Initial self test is "running".	Wait until the initial self test has completed.
Occlusion!	Tube or Veress needle has a blockage.	Localize the cause and open/eliminate the occlusion.
Leave obturator in place.	AirSeal initiating.	Leave obturator in place.

Error message	Cause	Remedy
Venting valve error! Restart the device. If the error occurs again call service!	Valve malfunction.	Restart the device. If the error occurs again, call service!
Electronic defective! Restart the device. If the error occurs again, call service!	Malfunctions in the electronic system.	Restart the device. If the error occurs again, call service!
Sensor error! Restart the device. If the error occurs again call service!	Sensor malfunction.	Restart the device. If the error occurs again, call service!
Device temperature error! Do not use device. Call service!	Temperature in device too high.	Serious defect. Do not use device. Call service!

Calibration error! The device must be re-calibrated. Call service!	The device is not calibrated properly.	The device must be re-calibrated. Call service!
Contamination! Fluid has penetrated the Control Unit internally! Ongoing surgery can be finished! Do not use device for next procedure! Contact technical service!	The device is contaminated with fluid.	The device must be checked by an authorized service technician or clearly marked with a label referring to the contamination and then twice enclosed in a safety foil, sealed, and returned to the manufacturer for repairs.
Venting valve error! Restart the device. If the error occurs again call service!	Valve malfunction.	Serious defect. Do not use device. Call service!

Critical message	Cause	Remedy
Check connection of Tri-Lumen Filtered Tube Set!	Tube set improperly connected.	Connect the tube set correctly.
Cannula occlusion! Ensure black lines on distal cannula are visible inside cavity!	Tube or cannula has a blockage.	Determine the cause for the occlusion. Check the positioning of the cannula.
Filter contaminated! Insert new tube set!	High fluid level of the fluid trap in the filter housing	Please insert a new tube set.
Filter contaminated! AirSeal insufflation stopped. Insert obturator into AirSeal cannula to maintain pressure! Change tube set to continue with AirSeal function!	High fluid level of the fluid trap in the filter housing during AirSeal mode.	Finish procedure and call service.
Filter contaminated! Smoke Evacuation function stopped! Change tube set to continue Smoke Evacuation function!	High fluid level of the fluid trap in the filter housing during Smoke Evacuation mode.	Finish procedure and call service.
Change gas bottle! Insufflation stopped!	Low insufflation gas supply (bottle or house supply) during insufflation.	Check gas supply.
Change gas bottle!	Low insufflation gas supply (bottle or house supply).	Check gas supply.
Change gas bottle! AirSeal Insufflation shut down in 100 s! Pressure limited to 12 mmHg!	Low insufflation gas supply (bottle or house supply) during AirSeal® insufflation.	Check gas supply.
Insufficient gas! Check gas supply and restart device!	Device has no gas supply.	Check gas supply.
Check gas supply!	Device has no gas supply.	Check gas supply.
Check gas supply! Insufflation stopped!	Low insufflation gas supply (bottle or house supply) during insufflation.	Check gas supply.
Check gas supply! AirSeal Insufflation will stop in 100 s! Pressure limited to 12 mmHg!	Low insufflation gas supply (bottle or house supply) during AirSeal® insufflation.	Check gas supply.
Gas level low. Prepare to change gas bottle.	Low insufflation gas supply (bottle or house supply) during insufflation.	Check gas supply.
Overpressure!	The actual pressure is 3 mmHg above the nominal pressure for longer than 5 seconds.	Determine the cause for exceeding the nominal pressure. Check the electronic displays of the device if overpressure exists for a longer period of time.
Overpressure! Venting active!	The actual pressure is above the nominal pressure.	Determine the cause for exceeding the nominal pressure. Check the electronic displays of the device if overpressure exists for a longer period of time.

13 Technical Data

Model		AirSeal® i.F.S. Intelligent Flow System		
Mains voltage range ¹ [V]		100V/115V/230V		
Supply frequency range [Hz]		EU: 230 V 50 Hz / 230 V 60 Hz US: 115 V 60 Hz / 100 V 60 Hz/ 100 V 50 Hz		
Fuse designation ²		2xT 10 A, 250 VAC, UL-recognized		
Power consumption		Current [A]	Voltage [V]	Power consumption [VA]
EU: 230 V 50 Hz				
	Normal operation	2.5	230	575
	Full load	4.0	230	920
EU: 230 V 60 Hz				
	Normal operation	2.9	230	667
	Full load	4.5	230	1035
US: 115 V 60 Hz				
	Normal operation	5.5	115	633
	Full load	7.6	115	874
US: 100 V 50 Hz				
	Normal operation	5.9	100	590
	Full load	6.9	100	690
US: 100 V 60 Hz				
	Normal operation	6.6	100	660
	Full load	7.8	100	780
Protection class (I, II, III)		I		
Application part of the type (B, BF, CF)		BF		
Protection type (IP code ³)		21		
Classification ⁴ (I, IIa, IIb, III)		IIa		
Conformity with the following standards:		DIN EN 60601-1:2007 / IEC 60601-1:2005 EN 60601-1-2:2007 / IEC 60601-1-2:2007		
Operating conditions		10 to 30 °C / 50 to 86 °F 30 to 75 % rel. humidity 700 to 1060 hPa air pressure 3000 m max. altitude above sea level for device use		
Storage and transportation conditions		-25 to +70 °C / -13 to +158°F 10 to 95 % rel. humidity 700 to 1060 hPa air pressure		
Max. sound level:		71 db		
Inlet pressure range		3.4 - 80 bar/49 - 1160 psi		
Flow volume (range)		40 sl/min.		
Adjustable values	Pressure [mmHg]	5-20		
	Flow [l/min]	1-40		
Accuracy ⁵	Pressure [mmHg]	±1 mmHg		
	Flow [l/min]	±2.5 sl/min.		
Accuracy ⁶	Pressure [mmHg]	±1 mmHg		
	Flow [l/min]	±2.5 sl/min.		

Dimensions	Width x Height x Depth [mm ³]	420 x 220 x 470
Weight [kg]		26

Interfaces:	Service interface (USB port)
	Mains power connection (IEC-60320-1 C14)
Software version:	See User Menu

1. Identification according to IEC 61293; example: 100-240V or 110V/220V
2. Identification necessary if fuse holder is a contact part; information about voltage, current, trigger speed and cutoff capacity required; e.g.: 2x T 3,15 AH, 250 V, UL-recognized
3. See DIN EN 60529; part 9
4. Acc. to 93/42/EEC
5. Formerly "repetition accuracy" - describes the reproducibility of a measurement ("internal precision")
6. Degree of match between displayed (actual) and correct (nominal) value; typical info "x% Full Scale"

14 AirSeal®iFS Warranty

Standard Warranty

SurgiQuest, Inc. (the "Company") warrants this product against defects in materials and workmanship to the registered owner at the time of purchase. All components are covered by the warranty for a period of one year from the date of purchase. The Company will repair such defects in materials and workmanship during the warranty period at no charge to the customer.

This warranty does not apply to any unit that has been subject to misuse, neglect, improper installation or that has been altered, adjusted, or tampered with by any person other than the Company's authorized personnel.

If upon examination by authorized service personnel, it is determined that a malfunction is due to misuse, abuse or any other circumstance set forth in the preceding paragraph, the warranty provisions will not apply. In such a case, an estimate of the cost of repair work will be provided to the customer prior to servicing and repairing the unit.

During the warranty period for covered repairs, SurgiQuest, Inc. shall bear the cost of shipping for both the defective unit and any replacement equipment provided. After the warranty period or with respect to repairs not covered by this warranty, however, the customer is responsible for returning the defective equipment to the company authorized repair location at its own expense. The Company or its representative will service the unit, repair or replace any defective parts thereof, and return the unit.

If, upon examination, it is determined that the fault is not covered by this warranty, the repairs will be billed to the customer as out-of-warranty repairs. Devices repaired under the Company's standard repair program will remain covered by this warranty for the longer of thirty days or the balance of the original product warranty.

The warranty and the remedies as set forth herein are exclusive and in lieu of all other warranties, remedies, obligations and liabilities of SurgiQuest, Inc., expressed or implied, including the implied warranties of merchantability and fitness for a particular purpose. In no event shall SurgiQuest, Inc. be liable for any breach of warranty in any amount exceeding the purchase price of the product.

No agent, employee or representative of SurgiQuest, Inc. has the authority to bind the Company to any other warranty, affirmation, or representation concerning this device.

This warranty is valid only to the original purchaser of SurgiQuest, Inc. products directly from SurgiQuest, Inc. or from a SurgiQuest, Inc. authorized agent. The warranty cannot be transferred or assigned by the original purchaser.

Extended Warranty

Purchase of the "Extended Warranty" extends the "Standard Warranty" an additional twelve months for a total protection period of two years. Purchase must be completed within the "Standard Warranty" coverage period.

Premium Warranty

Purchase of the "Premium Warranty" extends the "Standard Warranty" an additional twenty-four months for a total protection period of three years. Purchase must be completed within the "Standard Warranty" coverage period.

Contact your local SurgiQuest representative for additional information on the Extended and Premium Warranties

Service and Claims:**EN**

If service is needed either during or after the warranty period:

- In the United States, contact SurgiQuest, Inc.
at +1-414-434-6634 or call your local SurgiQuest, Inc. sales representative.
- Package all the components carefully in the original shipping container. If the original shipping container is not available, an equivalent shipping container must be obtained from SurgiQuest.
- Print the Return Authorization supplied by SurgiQuest, Inc. and prepare for shipment pickup at scheduled time and location.

15 Test Log

15.1 Test Log

[illegible]

Please fill out and include this form when returning the device:

--

Street:

--

[illegible]

Date _____

Index**A**

Access to the service menu 29

Authorized trained personnel 32, 32

C

Care and maintenance 4

Certification 32

Contamination 4, 27

Contraindications 7

D

Delivery inspection 11

E

ESD (Electrostatic Discharge) precautionary measures 37

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Copyright © SURGIQUEST, Inc.